



**syrex**

**Disposable Syringe**

ISO13485:2016 **10ml Luer Slip** CE 0123

Manufactured for **Balmer Healthcare Limited**  
by Jiangsu Kanghua Medical Equipment Co., Ltd.  
Address: Sanhekou, Changzhou, Jiangsu, China. P.C.: 213115  
Tel: 0086 519 88676715 88670071  
Fax: 0086 519 88676688  
E-mail: kanghua@chinakang-hua.com  
Http: //www.kanghuamedical.com



STERILE EO

- Luer Slip Syringe 2-part with needle, latex free
- Non-bacteria, non-toxic, non-pyrogenic;
- Don't use if the package is open or damaged;



27L 1X SN: 1845  
U 100 216  
10031  
19031  
2019031  
2019031  
2019031

DISPOSABLE SYRINGE **syrex**



1  
2  
3  
4  
5  
6  
7  
8  
9  
10  
11  
12  
cc/mL





100  
PCS

1. Sterility not guaranteed if the individual package is punctured, opened or Damaged.
2. Do not use if the individual package is stained, dirty or damaged.
3. Single use only.
4. Discard the syringe properly after use.
5. Store in a clean and dry place at a temperature between 15°C and 25°C. Store away from sunlight.

# Disposable Syringe

INSULIN (U-100)

CE 0123

0.5ml 29G x 1/2" STERILE EO ISO13485:2016


**Manufactured for Balmer Healthcare Limited**  
 by Jiangsu Kanghua Medical Equipment Co., Ltd.  
 Address: Sanhekou, Changzhou, Jiangsu, China. P.C.: 213115  
 Tel: 0086 519 88676715 88670071 Fax: 0086 519 88678968  
 E-mail: kanghua@chinakang-hua.com www.kanghuamedical.com

# syrex®




# syrex®

**Disposable Syringe**  
 INSULIN (U-100)  
 0.5ml 29G x 1/2"  
 100 PCS



# syrex®



## IV Set with Flow Regulator ( Adult )

# Infurex

### TO USE INFUSION SET

- Close the roller clamp
- Remove the cover from plastic lancet
- Punch the lancet into the fluid container
- Open the roller clamp and expel the air bubbles
- Close the roller clamp
- Adjust the flow rate
- 1ml = 20 drops
- Luer Lock Connector

- Do not store at extreme temperature and humidity
- NON-TOXIC
- NON-PYROGENIC
- Discard after single use
- Sterilized by EO



STERILE EO



Please use this pouch to dispose the enclosed medical device after use

LOT



Manufactured for Prospectus Therapeutics Ltd, P.O. Box 54000-00200, Email: [sales@prospectus.co.ke](mailto:sales@prospectus.co.ke),  
Website: [www.prospectus.co.ke](http://www.prospectus.co.ke) by Hospital & Homecare IMP & EXP. Co. Ltd Nanjing





# Lightproof IV Set with Flow Regulator ( Adult )

## Infurex

### TO USE INFUSION SET

- Close the roller clamp
- Remove the cover from plastic lancet
- Punch the lancet into the fluid container
- Open the roller clamp and expel the air bubbles
- Close the roller clamp
- Adjust the flow rate
- 1ml = 20 drops
- Luer Lock Connector

- Do not store at extreme temperature and humidity
- NON-TOXIC
- NON-PYROGENIC
- Discard after single use
- Sterilized by EO



STERILE EO



Please use this pouch to dispose the enclosed medical device after use

LOT



Manufactured for Prospectus Therapeutics Ltd, P.O. Box 54000-00200, Email: sales@prospectus.co.ke,  
Website: www.prospectus.co.ke by Hospital & Homecare (MP & EXP. Co. Ltd Naringog

